


ABB Stock Lenses



Stock Lens Portfolio

Easy Reference Guide

CR-39 Plastic

Mfr/ Vendor	WDS Code	Description	Blank Size mm	Power Availability	Scratch Resistant	Anti- Reflective (AR)	Hydrophobic Water Resistant	Oleophobic Grease Resistant	Anti- Static
ABB	AC3	ABB PRIME AR	65/70	Sph: +6.00 to -6.00 Cyl: +6.00 to -6.00 to 3.00 cyl	✓	✓	✓		
SOMO	UC39	UltraClear AR with ITO	65 (plus powers) 70 (minus powers)	Sph: +4.00 to -6.00 Cyl: +4.00 to -6.00 to 2.00 cyl	✓	✓	✓	✓	✓
Essilor	CS39	Crizal Sapphire HR	70	Sph: +4.00 to -6.00 Cyl: +4.00 to -4.00 to 2.00 cyl, -4.25 to -5.75 Combined power cannot exceed -6.00	✓	✓	✓	✓	✓
Shamir	GLX39	Glacier Expression AR <i>NEW</i>	65 (plus powers) 70 (minus powers)	Sph: +3.00 to -4.00 Cyl: +3.00 to -3.00 to 2.00 cyl, -3.25 to 1.75 cyl, -3.50 to 1.50 cyl, -3.75 to 1.25 cyl, -4.00 to -1.00 cyl	✓	✓	✓	✓	✓
Essilor	CRO	Crizal Rock	70	Sph: +4.00 to -6.00 Cyl: +4.00 to -4.00 to 2.00 cyl, -4.25 to -5.75 Combined power cannot exceed -6.00	✓	✓	✓	✓	✓
Shamir	GLP39	Glacier Plus AR <i>NEW</i>	65 (plus powers) 70 (minus powers)	Sph: +3.00 to -4.00 Cyl: +3.00 to -4.00 to 2.00 cyl	✓	✓	✓	✓	✓
Essilor	CEO	Crizal Easy Pro	70	Sph: +4.00 to -6.00 Cyl: +4.00 to -4.00 to 2.00 cyl, -4.25 to -5.75 Combined power cannot exceed -6.00	✓	✓	✓	✓	✓
SOMO	1539	SOMO AR	65 (plus powers) 70 (minus powers)	Sph: +4.00 to -6.00 Cyl: +4.00 to -6.00 to 2.00 cyl	✓	✓	✓		
SOMO	HD39	SOMO Hard Coat <i>TINTABLE</i>	65 (plus powers) 70 (minus powers)	Sph: +4.00 to -6.00 Cyl: +4.00 to -6.00 to 2.00 cyl	✓				
Essilor	E70C	Essilor CR-39 Flat True Tint 70mm <i>TINTABLE</i>	70	Sph: +6.00 to -8.00 Cyl: +6.00 to -5.00 to 3.00 cyl, -5.25 to -7.75 Combined power cannot exceed -8.00	✓				
Essilor	E75C	Essilor CR-39 Flat True Tint 75mm <i>TINTABLE</i>	75	Sph: +4.00 to -6.00 Cyl: +4.00 to -4.00 to 2.00 cyl, -4.25 to -5.75 Combined power cannot exceed -6.00	✓				
SOMO	HDUC	Uncoated 65/70mm <i>TINTABLE</i>	65 (plus powers) 70 (minus powers)	Sph: +4.00 to -6.00 Cyl: +4.00 to -6.00 to 2.00 cyl					

Specialty

Mfr/ Vendor	WDS Code	Description	Blank Size mm	Power Availability	Scratch Resistant	Anti- Reflective (AR)	Hydrophobic Water Resistant	Oleophobic Grease Resistant	Anti- Static
SOMO	BOOST	Hue ColorBoost Polycarbonate AR 2pk <i>Color Option: Green Indoor, Green Outdoor</i>	78	Plano Only <i>Base Curves Available: 4, 6, and 8</i>	✓	✓			
Shore	PAFH	Polycarbonate Hard Coat Anti-fog	65 (plus powers) 70 (minus powers)	Sph: +4.00 to -6.00 Cyl: +4.00 to -6.00 to 2.00 cyl	✓				
Shore	SHFT	Flat Top 28 CR-39 <i>Right/Left</i>	70	Dist: 0.00 to +3.00, Add Powers +1.00 to +3.00 Dist: -0.25 to -2.00, Add Powers +1.00 to +2.00					
Conant	C39	CR-39 Clear & Tint 25 pair pack <i>Color Option: Clear, Gray, G15, Brown</i>	74	Plano Only					

Polycarbonate									
Mfr/ Vendor	WDS Code	Description	Blank Size mm	Power Availability	Scratch Resistant	Anti- Reflective (AR)	Hydrophobic Water Resistant	Oleophobic Grease Resistant	Anti- Static
ABB	ACP	ABB PRIME AR	65/70	Sph: +6.00 to -6.00 Cyl: +6.00 to -6.00 to 2.00 cyl	✓	✓	✓		
SOMO	UCPO	UltraClear AR with ITO	68 (plus powers) 72 (minus powers)	Sph: +4.00 to -8.00 Cyl: +4.00 to -8.00 to 2.00 cyl (-6.25 to -8.00 are Aspheric)	✓	✓	✓	✓	✓
Essilor	CSP	Crizal Sapphire HR	65/70	Sph: +4.00 to -8.00 Cyl: +4.00 to -6.00 to 2.00 cyl; -6.25 to -8.00 Combined power cannot exceed -8.00	✓	✓	✓	✓	✓
Shamir	GLXPY	Glacier Expression AR <i>NEW</i>	65/70	Sph: +4.00 to -6.00 Cyl: +4.00 to -6.00 to 2.00 cyl Combined power cannot exceed -6.00	✓	✓	✓	✓	✓
Essilor	CRP	Crizal Rock	65/70	Sph: +4.00 to -8.00 Cyl: +4.00 to -6.00 to 2.00 cyl; -6.25 to -8.00 Combined power cannot exceed -8.00	✓	✓	✓	✓	✓
Shamir	GLPPY	Glacier Plus AR <i>NEW</i>	65/70	Sph: +4.00 to -6.00 Cyl: +4.00 to -6.00 to 2.00 cyl Combined power cannot exceed -6.00	✓	✓	✓	✓	✓
Essilor	CEP	Crizal Easy Pro	65/70	Sph: +4.00 to -8.00 Cyl: +4.00 to -6.00 to 2.00 cyl; -6.25 to -8.00 Combined power cannot exceed -8.00	✓	✓	✓	✓	✓
SOMO	15PO	SOMO AR	68 (plus powers) 72 (minus powers)	Sph: +4.00 to -8.00 Cyl: +4.00 to -8.00 to 2.00 cyl (-6.25 to -8.00 are Aspheric)	✓	✓	✓		
SOMO	HDPO	SOMO Hard Coat	68 (plus powers) 72 (minus powers)	Sph: +4.00 to -8.00 Cyl: +4.00 to -8.00 to 2.00 cyl (-6.25 to -8.00 are Aspheric)	✓				
SOMO	STPO	Super Tint Hard Coat <i>TINTABLE</i>	68 (plus powers) 72 (minus powers)	Sph: +4.00 to -6.00 Cyl: +4.00 to -6.00 to 2.00 cyl	✓				

Polycarbonate Aspheric									
Mfr/ Vendor	WDS Code	Description	Blank Size mm	Power Availability	Scratch Resistant	Anti- Reflective (AR)	Hydrophobic Water Resistant	Oleophobic Grease Resistant	Anti- Static
SOMO	UCPA	UltraClear AR with ITO <i>w/Nano No-Slip</i>	60 (+4.25 to +6.00) 68 (+0.25 to +4.00) 72 (0.00 to -8.00)	Sph: +6.00 to -8.00 Cyl: +0.25 to +3.00 to 3.00 cyl, +3.25 to +6.00 to 2.00 cyl, 0.00 to -6.00 to 3.00 cyl, -6.25 to -8.00 to 2.00 cyl	✓	✓	✓	✓	✓
Essilor	CRPA	Crizal Rock	65/70	Sph: +4.00 to -8.00 Cyl: +4.00 to -6.00 to 2.00 cyl; -6.25 to -8.00 Combined power cannot exceed -8.00	✓	✓	✓	✓	✓
Essilor	EAAR	Essilor AR	65 (plus powers) 70 (minus powers; cyl to -2.00) 65 (minus powers; cyl -2.25 to -4.00)	Sph: +6.00 to -10.00 Cyl: +6.00 to -6.00 to 4.00 cyl, -6.25 to -9.75 Combined power cannot exceed -10.00	✓	✓	✓		
Essilor	EAA	Essilor Hard Coat <i>TINTABLE</i>	65 (plus powers) 70 (minus powers; cyl to -2.00) 65 (minus powers; cyl -2.25 to -4.00)	Sph: +6.00 to -10.00 Cyl: +6.00 to -6.00 to 4.00 cyl, -6.25 to -9.75 Combined power cannot exceed -10.00	✓				

Trivex									
Mfr/ Vendor	WDS Code	Description	Blank Size mm	Power Availability	Scratch Resistant	Anti- Reflective (AR)	Hydrophobic Water Resistant	Oleophobic Grease Resistant	Anti- Static
Shamir	GLPTX	Glacier Plus AR <i>NEW</i>	65 (plus powers) 70 (minus powers)	Sph: +4.00 to -6.00 Cyl: +4.00 to -6.00 to 2.00 cyl	✓	✓	✓	✓	✓
TOG	TAR	Trivex AR	65/70	Sph: +3.00 to -6.00 Cyl: +3.00 to -6.00 to 2.00 cyl	✓	✓	✓	✓	✓
Seiko/Hoya	HTXA	Hoya AR	65 (plus powers) 70 (minus powers)	Sph: +3.00 to -8.00 Cyl: +3.00 to -7.00 to 2.00 cyl, -7.25 to 1.50 cyl, -7.50 to 1.00 cyl, -7.75 to 0.50 cyl	✓	✓	✓	✓	✓
Seiko/Hoya	HTX	Hoya Hard Coat	65 (plus powers) 70 (minus powers)	Sph: +4.00 to -6.00 Cyl: +4.00 to -6.00 to 2.00 cyl	✓				

1.60 Spherical									
Mfr/ Vendor	WDS Code	Description	Blank Size mm	Power Availability	Scratch Resistant	Anti- Reflective (AR)	Hydrophobic Water Resistant	Oleophobic Grease Resistant	Anti- Static
SOMO	S6SA	SOMO AR	65 (+2.25 to +6.00) 70 (+0.25 to +2.00) 75 (pl to -6.00) 70 (-6.25 to -8.00)	Sph: +6.00 to -8.00 Cyl: +0.25 to +6.00 to 2.00 cyl, 0.00 to -6.00 to 4.00 cyl, -6.25 to -8.00 to 2.00 cyl	✓	✓	✓		
SOMO	S6S	SOMO Hard Coat	65 (+2.25 to +6.00) 70 (+0.25 to +2.00) 75 (pl to -6.00) 70 (-6.25 to -8.00)	Sph: +6.00 to -8.00 Cyl: +0.25 to +6.00 to 2.00 cyl, 0.00 to -6.00 to 4.00 cyl, -6.25 to -8.00 to 2.00 cyl	✓				

1.67 Aspheric

Mfr/ Vendor	WDS Code	Description	Blank Size mm	Power Availability	Scratch Resistant	Anti- Reflective (AR)	Hydrophobic Water Resistant	Oleophobic Grease Resistant	Anti- Static
ABB	AC6	ABB PRIME AR	65/70/75	Sph: +6.00 to -15.00 Cyl: +6.00 to -12.00 to 2.00 cyl -12.50 to -15.00 (0.50D) plano cyl	✓	✓	✓		
SOMO	UC7A	UltraClear AR with ITO	75 (+0.25 to +2.00); 70 (+2.25 to +4.00) 65 (+4.25 to +6.00); 75 (pl to -6.00) 70 (-6.25 to -10.00)	Sph: +6.00 to -10.00 Cyl: +0.25 to +6.00 to 2.00 cyl, 0.00 to -8.00 to 3.00 cyl, -8.25 to -10.00 to 2.00 cyl	✓	✓	✓	✓	✓
Essilor	CS67	Crizal Sapphire HR	65 - 70 - 75 - 80	Sph: -3.00 to -12.00 Cyl: -3.00 to -10.25 to 2.00 cyl, -10.50 to -12.00 plano cyl	✓	✓	✓	✓	✓
Shamir	GLX67	Glacier Expression AR <i>NEW</i>	65 (plus powers) 80 (minus powers)	Sph: +4.00 to -6.00 Cyl: +4.00 to -6.00 to 2.00 cyl	✓	✓	✓	✓	✓
Essilor	CR67	Crizal Rock	65 - 70 - 75 - 80*	Sph: +6.00 to -12.00 Cyl: +6.00 to -10.25 to 2.00 cyl, -10.50 to -12.00 plano cyl	✓	✓	✓	✓	✓
Shamir	GLP67	Glacier Plus AR <i>NEW</i>	65 (plus powers) 75 (minus powers)	Sph: +4.00 to -6.00 Cyl: +4.00 to -6.00 to 2.00 cyl	✓	✓	✓	✓	✓
Essilor	CE67	Crizal Easy Pro	65 - 70 - 75 - 80*	Sph: +6.00 to -12.00 Cyl: +6.00 to -10.00 to 2.00 cyl, -10.25 to -12.00 plano cyl	✓	✓	✓	✓	✓
SOMO	HD7A	SOMO Hard Coat	75 (+0.25 to +2.00) 75 (pl to -6.00) 70 (-6.25 to -10.00)	Sph: +2.00 to -10.00 Cyl: +0.25 to +2.00 to 2.00 cyl, 0.00 to -8.00 to 3.00 cyl, -8.25 to -10.00 to 2.00 cyl	✓				

1.74 Aspheric

Mfr/ Vendor	WDS Code	Description	Blank Size mm	Power Availability	Scratch Resistant	Anti- Reflective (AR)	Hydrophobic Water Resistant	Oleophobic Grease Resistant	Anti- Static
SOMO	UC74	UltraClear AR with ITO	75 (-2.25 to -8.00) 70 (-8.25 to -12.00) 65 (-12.50 to -15.00)	Sph: -2.25 to -15.00 (-10.50 to -15.00 0.50D only) Cyl: -2.25 to -8.00 to 4.00 cyl, -8.25 to -10.00 to 3.00 cyl, -10.50 to -12.00 to 2.00 cyl (0.50D only)	✓	✓	✓	✓	✓
Essilor	CS74	Crizal Sapphire HR	70/75	Sph: -3.00 to -12.00 Cyl: +6.00 to -10.25 to 2.00 cyl, -10.50 to -12.00 plano cyl	✓	✓	✓	✓	✓

Photochromic CR-39

Mfr/ Vendor	WDS Code	Description	Blank Size mm	Power Availability	Scratch Resistant	Anti- Reflective (AR)	Hydrophobic Water Resistant	Oleophobic Grease Resistant	Anti- Static
Hoya	S2CA	Sensity 2 EX3+ AR <i>Color Option: Gray</i>	65/70	Sph: +3.00 to -4.00 Cyl: +3.00 to -4.00 to 2.00 cyl	✓	✓	✓	✓	✓
Essilor	TGS3S	Crizal Sapphire HR - Transitions GEN S <i>Color Option: Gray</i>	65/70	Sph: +4.00 to -4.00 Cyl: +4.00 to -2.00 to 2.00 cyl; -2.25 to -4.00 Combined power cannot exceed -4.00	✓	✓	✓	✓	✓
Essilor	TGS3R	Crizal Rock - Transitions GEN S <i>Color Option: Gray</i>	65/70	Sph: +4.00 to -4.00 Cyl: +4.00 to -2.00 to 2.00 cyl; -2.25 to -4.00 Combined power cannot exceed -4.00	✓	✓	✓	✓	✓
Essilor	TGS3E	Crizal Easy Pro - Transitions GEN S <i>Color Option: Gray</i>	65/70	Sph: +4.00 to -4.00 Cyl: +4.00 to -2.00 to 2.00 cyl; -2.25 to -4.00 Combined power cannot exceed -4.00	✓	✓	✓	✓	✓
Hoya	S2CH	Sensity 2 Hard Coat <i>Color Option: Gray</i>	65/70	Sph: +3.00 to -4.00 Cyl: +3.00 to -4.00 to 2.00 cyl	✓				
Essilor	TGS39	Transitions GEN S Hard Coat <i>Color Options: Gray, Brown</i>	65/70	Sph: +4.00 to -4.00 Cyl: +4.00 to -2.00 to 2.00 cyl; -2.25 to -4.00 Combined power cannot exceed -4.00	✓				

Photochromic Polycarbonate

Mfr/ Vendor	WDS Code	Description	Blank Size mm	Power Availability	Scratch Resistant	Anti- Reflective (AR)	Hydrophobic Water Resistant	Oleophobic Grease Resistant	Anti- Static
Hoya	S2PA	Sensity 2 EX3+ AR <i>Color Option: Gray</i>	65/70	Sph: +3.00 to -6.00 Cyl: +3.00 to -6.00 to 2.00 cyl	✓	✓	✓	✓	✓
Essilor	TGSPS	Crizal Sapphire HR - Transitions GEN S <i>Color Option: Gray</i>	65/70	Sph: +4.00 to -6.00 Cyl: +4.00 to -4.00 to 2.00 cyl; -4.25 to -6.00 Combined power cannot exceed -6.00	✓	✓	✓	✓	✓
Essilor	TGSPR	Crizal Rock - Transitions GEN S <i>Color Option: Gray</i>	65/70	Sph: +4.00 to -6.00 Cyl: +4.00 to -4.00 to 2.00 cyl; -4.25 to -6.00 Combined power cannot exceed -6.00	✓	✓	✓	✓	✓
Shamir	GLPPT	Glacier Plus AR - Transitions GEN S <i>Color Option: Gray</i> NEW	65 (plus powers) 70 (minus powers)	Sph: +4.00 to -6.00 Cyl: +4.00 to -6.00 to 2.00 cyl Combined power cannot exceed -6.00	✓	✓	✓	✓	✓
Essilor	TGSPE	Crizal Easy Pro - Transitions GEN S <i>Color Option: Gray</i>	65/70	Sph: +4.00 to -6.00 Cyl: +4.00 to -4.00 to 2.00 cyl; -4.25 to -6.00 Combined power cannot exceed -6.00	✓	✓	✓	✓	✓
VisionEase	CPPAAR	Choice Photo AR Aspheric <i>Color Option: Gray and Brown</i> NEW	65 (plus powers) 75 (minus powers)	<i>Color Option: Gray</i> Sph: +2.00 to -4.00 Cyl: +2.00 to -4.00 to 2.00 cyl <i>Color Option: Brown</i> Sph: +2.00 to -4.00 Cyl: +2.00 to -4.00 +1.25 to +2.00 & -3.25 to -4.00 sph only +0.25 to +1.00 & -1.25 to -3.00 cyl to 1.50 0.00 to -1.00 cyl to -2.00	✓	✓			
Hoya	S2PH	Sensity 2 Hard Coat <i>Color Option: Gray</i>	65 (+1.25 to +4.00) 70 (+1.00 to -6.00)	Sph: +4.00 to -6.00 Cyl: +4.00 to -6.00 to 2.00 cyl	✓				
Essilor	TGSPY	Transitions GEN S Hard Coat <i>Color Options: Gray, Brown</i>	65/70	Sph: +4.00 to -6.00 Cyl: +4.00 to -4.00 to 2.00 cyl; -4.25 to -6.00 Combined power cannot exceed -6.00	✓				
VisionEase	CPPAHC	Choice Photo Hard Coat Aspheric <i>Color Option: Gray and Brown</i> NEW	65 (plus powers) 75 (minus powers)	<i>Color Option: Gray</i> Sph: +2.00 to -4.00 Cyl: +2.00 to -4.00 to 2.00 cyl <i>Color Option: Brown</i> Sph: +2.00 to -4.00 Cyl: +2.00 to -4.00 +1.25 to +2.00 & -3.25 to -4.00 sph only +0.25 to +1.00 & -1.25 to -3.00 cyl to 1.50 0.00 to -1.00 cyl to -2.00	✓				

Photochromic Trivex									
Mfr/ Vendor	WDS Code	Description	Blank Size mm	Power Availability	Scratch Resistant	Anti- Reflective (AR)	Hydrophobic Water Resistant	Oleophobic Grease Resistant	Anti- Static
Hoya	S2TA	Sensity 2 EX3+ AR <i>Color Option: Gray</i>	65 (plus powers) 70 (minus powers)	Sph: +2.00 to -6.00 Cyl: +2.00 to -4.00 to 2.00 cyl, -4.25 to -6.00 Combined power cannot exceed -6.00	✓	✓	✓	✓	✓
Hoya	S2TH	Sensity 2 Hard Coat <i>Color Option: Gray</i>	65 (plus powers) 70 (minus powers)	Sph: +2.00 to -6.00 Cyl: +2.00 to -4.00 to 2.00 cyl, -4.25 to -6.00 Combined power cannot exceed -6.00	✓				

Blue Blocking AR CR39									
Mfr/ Vendor	WDS Code	Description	Blank Size mm	Power Availability	Scratch Resistant	Anti- Reflective (AR)	Hydrophobic Water Resistant	Oleophobic Grease Resistant	Anti- Static
ABB	AC3B	ABB PRIME BLUE AR (UV 420 in-mass)	65/70	Sph: +6.00 to -6.00 Cyl: +6.00 to -6.00 to 2.00 cyl	✓	✓	✓		
Essilor	CPOF	Previncia Crizal AR Orma Flat (CR-39)	70	Sph: +4.00 to -6.00 Cyl: +4.00 to -4.00 to 2.00 cyl, -4.25 to -5.75 Combined power cannot exceed -6.00	✓	✓	✓	✓	✓

Blue Blocking AR Polycarbonate									
Mfr/ Vendor	WDS Code	Description	Blank Size mm	Power Availability	Scratch Resistant	Anti- Reflective (AR)	Hydrophobic Water Resistant	Oleophobic Grease Resistant	Anti- Static
ABB	ACPB	ABB PRIME BLUE AR (UV 420 in-mass)	65/70	Sph: +6.00 to -6.00 Cyl: +6.00 to -6.00 to 2.00 cyl	✓	✓	✓		
SOMO	UCPB	UltraClear AR with ITO Blue Block	68 (plus powers) 72 (minus powers)	Sph: +4.00 to -6.00 Cyl: +4.00 to -6.00 to 2.00 cyl	✓	✓	✓	✓	✓
Essilor	CPU	Previncia Crizal AR	65 (+1.00 to +4.00) 70 (+0.25 to +0.75) 70 (pl to -8.00)	Sph: +4.00 to -8.00 Cyl: +4.00 to +6.00 to 2.00 cyl, -6.25 to -7.75 Combined power cannot exceed -8.00	✓	✓	✓	✓	✓
Essilor	CSBC	Crizal Sapphire HR Blue Capture	65/70	Sph: +4.00 to -8.00 Cyl: +4.00 to +6.00 to 2.00 cyl Combined power cannot exceed -8.00	✓	✓	✓	✓	✓
Shamir	GLPPB	Glacier Plus Blue Zero AR <i>NEW</i>	65 (plus powers) 70 (minus powers)	Sph: +4.00 to -6.00 Cyl: +4.00 to -6.00 to 2.00 cyl	✓	✓	✓	✓	✓
Essilor	CRBC	Crizal Rock Blue Capture	65/70	Sph: +4.00 to -8.00 Cyl: +4.00 to +6.00 to 2.00 cyl Combined power cannot exceed -8.00	✓	✓	✓	✓	✓

Blue Blocking AR Trivex									
Mfr/ Vendor	WDS Code	Description	Blank Size mm	Power Availability	Scratch Resistant	Anti- Reflective (AR)	Hydrophobic Water Resistant	Oleophobic Grease Resistant	Anti- Static
Shamir	GLPTB	Glacier Plus Blue Zero AR <i>NEW</i>	65 (plus powers) 70 (minus powers)	Sph: +4.00 to -6.00 Cyl: +4.00 to -6.00 to 2.00 cyl	✓	✓	✓	✓	✓
TOG	TAB	Bluloc AR (UV 420 in-mass)	65/70	Sph: +3.00 to -6.00 Cyl: +3.00 to -6.00 to 2.00 cyl	✓	✓	✓	✓	✓
Seiko/Hoya	HTRA	Recharge AR	65 (plus powers) 70 (minus powers)	Sph: +2.00 to -6.00 Cyl: +2.00 to -6.00 to 2.00 cyl	✓	✓	✓	✓	✓

Blue Blocking AR Mid Index									
Mfr/ Vendor	WDS Code	Description	Blank Size mm	Power Availability	Scratch Resistant	Anti- Reflective (AR)	Hydrophobic Water Resistant	Oleophobic Grease Resistant	Anti- Static
SOMO	UC6B	UltraClear AR with ITO 1.60 Aspheric	65 (+0.25 to +6.00) 75 (PL to -6.00) 70 (-6.25 to -10.00)	Sph: +6.00 to -10.00 Cyl: +0.25 to +6.00 to 2.00 cyl, 0.00 to -8.00 to 3.00 cyl, -8.25 to -10.00 to 1.00 cyl	✓	✓	✓	✓	✓
EST	BTI	BluTech High Impact AR Indoor 1.56 6 Base <i>Color Option: Brown Tint</i>	70/75	Sph: +2.00 to -3.00 Cyl: +2.00 to -3.00 to 1.00 cyl	✓	✓	✓	✓	✓
EST	BTIP	BluTech High Impact AR Indoor 1.56 6 Base <i>Color Option: Brown Tint</i>	75	Plano Only	✓	✓	✓	✓	✓

Blue Blocking AR Hi Index									
Mfr/ Vendor	WDS Code	Description	Blank Size mm	Power Availability	Scratch Resistant	Anti- Reflective (AR)	Hydrophobic Water Resistant	Oleophobic Grease Resistant	Anti- Static
ABB	AC6B	ABB PRIME BLUE AR 1.67 Aspheric (UV 420 in-mass)	65/70/75	Sph: +6.00 to -15.00 Cyl: +6.00 to -12.00 to 2.00 cyl -12.50 to -15.00 (0.50D) plano cyl	✓	✓	✓		
SOMO	UB67	UltraClear AR with ITO 1.67 Aspheric (UV 420 in-mass)	65/70/75	Sph: +6.00 to -12.00 Cyl: Plano to -8.00 to 4.00 cyl +0.25 to +6.00 & -8.25 to -12.00 to 2.00 cyl	✓	✓	✓	✓	✓
SOMO	UB74	UltraClear AR with ITO 1.74 Aspheric (UV 420 in-mass)	70/75	Sph: -4.00 to -11.00 (0.50D after -10.00) Cyl: -4.00 to -10.00 to 2.50 cyl -10.50 to -11.00 to 2.00 cyl	✓	✓	✓	✓	✓

Polarized Plano									
Mfr/ Vendor	WDS Code	Description	Blank Size mm	Power Availability	Scratch Resistant	Anti- Reflective (AR)	Hydrophobic Water Resistant	Oleophobic Grease Resistant	Anti- Static
PLASTIC									
Younger	YPP	Plastic 6 Base <i>Color Options: Brown, Gray Green, Gray</i>	76	Plano Only					
POLYCARBONATE									
KBCO	KBPM	Polycarbonate Mirror Standard Colors 6 Base <i>Mirror Color Options w/ Brown Lens:</i> Gold (GB) <i>Mirror Color Options w/ Gray Lens:</i> Silver (SG), Blue (BG), Green (GG)	71	Plano Only	✓				
KBCO	KBPMP	Polycarbonate Mirror Premium Colors 6 Base <i>Mirror Color Options w/ Brown Lens:</i> Red (RB) <i>Mirror Color Options w/ Gray Lens:</i> Ice Blue (IG), Purple (RG), Pink (PG), Chrome (CH), Black (BK), Rose Gold (RS), Champagne (CP)	71	Plano Only	✓				
KBCO	KBP	Polycarbonate 6 Base <i>Color Options: Brown, Gray, True Gray</i>	71	Plano Only	✓				

COLOR CODES:

- BR = Brown
- GG = Gray Green
- GY = Gray
- TG = True Gray




ABB Contact Lens | ABB Labs | ABB Business Solutions

800.852.8089 | abboptical.com