















DELTA SCLERAL

Diagnostic Marks

Delta Scleral lenses have diagnostic marks to aid the fitter in identifying various fitting characteristics:

Code / Name	Prescription Designs	
	OD	OS
DS-DX Delta Scleral Dx	Flat Meridian Axis Lines Clear Single BLACK dot @apex Single BLACK dot @base Black Ink 	
DS Delta Scleral	Flat Meridian Axis Lines Clear Single BLACK dot @base Black Ink 	Flat Meridian Axis Lines Clear Double BLACK dots @base Black Ink 
DS-V Delta Scleral w/EV	Flat Meridian Axis Lines Clear Edge Vault inner perimeter (x3) Clear Single BLACK dot @base Black Ink 	Flat Meridian Axis Lines Clear Edge Vault inner perimeter (x3) Clear Double BLACK dots @base Black Ink 
DS-V2 Delta Scleral w/EVx2	Flat Meridian Axis Lines Clear Edge Vault inner perimeter (x6) Clear Single BLACK dot @base Black Ink 	Flat Meridian Axis Lines Clear Edge Vault inner perimeter (x6) Clear Double BLACK dots @base Black Ink 
DS-Q Delta Scleral w/QH	Quad Meridian Axis Lines Q1 Meridian with Line+Dot Clear Single BLACK dot @base Black Ink 	Quad Meridian Axis Lines Q1 Meridian with Line+Dot Clear Double BLACK dots @base Black Ink 
DS-QV Delta Scleral w/QH+EV	Quad Meridian Axis Lines Q1 Meridian with Line+Dot Clear Edge Vault inner perimeter (x3) Clear Single BLACK dot @base Black Ink 	Quad Meridian Axis Lines Q1 Meridian with Line+Dot Clear Edge Vault inner perimeter (x3) Clear Double BLACK dots @base Black Ink 
DS-QV2 Delta Scleral w/QH+EVx2	Quad Meridian Axis Lines Q1 Meridian with Line+Dot Clear Edge Vault inner perimeter (x6) Clear Single BLACK dot @base Black Ink 	Quad Meridian Axis Lines Q1 Meridian with Line+Dot Clear Edge Vault inner perimeter (x6) Clear Double BLACK dots @base Black Ink 