



DELTA SCLERAL

Extended Dx Reference Lenses

When the physical 18 lens Dx set won't reach far enough, Delta can still reach it.

15.5 Diameter	BC	Lens Sag @ 12.60mm	Lens Sag @ 15.50mm
19	13.00	1475	2600
20	11.50	1775	2900
21	9.80	2075	3200
22	9.00	2375	3500
1	8.50	2675	3800
2	8.10	2975	4100
3	7.80	3275	4400
4	7.50	3575	4700
5	7.20	3875	5000
6	6.95	4175	5300
23	6.75	4475	5600
24	6.55	4775	5900
25	6.35	5075	6200
26	6.20	5375	6500
27	6.05	5675	6800
16.5 Diameter	BC	Lens Sag @ 13.40	Lens Sag @ 16.50
28	13.50	1575	2900
29	12.00	1875	3200
30	10.00	2175	3500
31	9.10	2475	3800
7	8.70	2775	4100
8	8.30	3075	4400
9	8.00	3375	4700
10	7.65	3675	5000
11	7.40	3975	5300
12	7.10	4275	5600
32	6.90	4575	5900
33	6.70	4875	6200
34	6.55	5175	6500
35	6.35	5475	6800
36	6.20	5775	7100
17.5 Diameter	BC	Lens Sag @ 14.20	Lens Sag @ 17.50
37	14.00	1725	3200
38	12.20	2025	3500
39	10.80	2325	3800
40	9.60	2625	4100
13	9.00	2925	4400
14	8.55	3225	4700
15	8.15	3525	5000
16	7.80	3825	5300
17	7.50	4125	5600
18	7.25	4425	5900
41	7.00	4725	6200
42	6.80	5025	6500
43	6.60	5325	6800
44	6.45	5625	7100
45	6.30	5925	7400

Reference Lens # 1-18 are in the physical Dx Kit.
Empirically fit extra shallow or extra deep by starting with the Ref# and going from there.

