



DELTA SCLERAL

## Extended Dx Reference Lenses

When the physical 18 lens Dx set won't reach far enough, DELTA can still reach it.

15.5	19	20	21	22	1	2	3	4	5	6	23	24	25	26	27
BC	13.00	11.50	9.80	9.00	8.50	8.10	7.80	7.50	7.20	6.95	6.75	6.55	6.35	6.20	6.05
Sag@ 12.60	1475	1775	2075	2375	2675	2975	3275	3575	3875	4175	4475	4775	5075	5375	5675
Sag@ 15.50	2600	2900	3200	3500	3800	4100	4400	4700	5000	5300	5600	5900	6200	6500	6800

16.5	28	29	30	31	7	8	9	10	11	12	32	33	34	35	36
BC	13.50	12.00	10.00	9.10	8.70	8.30	8.00	7.65	7.40	7.10	6.90	6.70	6.55	6.35	6.20
Sag@ 13.40	1575	1875	2175	2475	2775	3075	3375	3675	3975	4275	4575	4875	5175	5475	5775
Sag@ 16.50	2900	3200	3500	3800	4100	4400	4700	5000	5300	5600	5900	6200	6500	6800	7100

17.5	37	38	39	40	13	14	15	16	17	18	41	42	43	44	45
BC	14.00	12.20	10.80	9.60	9.00	8.55	8.15	7.80	7.50	7.25	7.00	6.80	6.60	6.45	6.30
Sag@ 14.20	1725	2025	2325	2625	2925	3225	3525	3825	4125	4425	4725	5025	5325	5625	5925
Sag@ 17.50	3200	3500	3800	4100	4400	4700	5000	5300	5600	5900	6200	6500	6800	7100	7400



Reference Lens # 1-18 are in the physical Dx Kit.  
Empirically fit extra shallow or extra deep by starting with the Ref# and going from there.

