DELTA SCLERAL

ABB Optical Group Order Form


☐ Original Order ☐ Adjustment Order ☐ Reorder

Patient		
Account #		Ship to:
Fitter		
Date		

15.5			
Dx #	BC	Pwr	Sag
1	8.50	Plano	3.800
2	8.10	-1.00	4.100
3	7.80	-2.00	4.400
4	7.50	-3.00	4.700
5	7.20	-4.00	5.000
6	6.95	-5.00	5.300

16.5			
Dx #	BC	Pwr	Sag
7	8.70	Plano	4.100
8	8.30	-1.00	4.400
9	8.00	-2.00	4.700
10	7.65	-3.00	5.000
11	7.40	-4.00	5.300
12	7.10	-5.00	5.600

17.5			
Dx #	BC	Pwr	Sag
13	9.00	Plano	4.400
14	8.55	-1.00	4.700
15	8.15	-2.00	5.000
16	7.80	-3.00	5.300
17	7.50	-4.00	5.600
18	7.25	-5.00	5.900



1 and 2 are for original orders only

	Fitting Procedure	OD			OS		
1	**Select Diameter based on HVID**	11.4 < = 15.5	11.5 to 12.3 16.5	>= 12.4 17.5	11.4 <= 15.5	11.5 to 12.3 16.5	>= 12.4 17.5
2	**Select Dx Reference Lens #** Choose lens that has the best central clearance	1-6	7-12	13-18	1-6	7-12	13-18

Original order - specify deltas required compared to the Dx reference lens
Adjustment orders - specify additional deltas required compared to the previous Rx lens
Reorders - Specify additional deltas compared to previous finalized lens

3	**Central Clearance Delta** (+) for more central clearance (-) for less, in microns.							
4	**Modify Base Curve (optional)** Mid-peripheral clearance change, (+) for more, (-) for less, in microns Note: BC changes will NOT affect lens sag or central clearance.							
5	**Steep Sag Delta** (+) for more SAG @ landing in steep meridian (-) for less, in microns							
6	**Limbal Clearance Deltas – Flat & Steep** (+) for more limbal clearance (-) for less, in microns.	Flat/Shallow		Steep/Deep		Flat/Shallow		Steep/Deep	
7	**Haptic Sag Deltas – Flat & Steep** (+) for more Haptic Sag (deeper, tighter) (-) for less Haptic Sag (shallower, looser), in microns	Flat		Steep		Flat		Steep	
(or)	**Haptic Sag Deltas - Quad Haptic**	Q1-0	Q2-90	Q3-180	Q4-270	Q1-0	Q2-90	Q3-180	Q4-270
8	**Flat Meridian Axis** Note position of Dx lines (flat haptic)	ccw(+) cw(-)				ccw(+) cw(-)			
9	**Over Refraction**								
10	**Edge Vault 1**	dec	axis	diam	height	dec	axis	diam	height
	Edge Vault 2								
11	**Material / Color**	☐ Optimum Extra ☐ Optimum Extreme ☐ Optimum Infinite ☐ Hexa 100 ☐ Acuity 100 ☐ Acuity 200 ☐ Boston XO				**OD** ☐ Clear ☐ Ice Blue		**OS** ☐ Clear ☐ Ice Blue	
12	**Hydrapeg**	☐ Hydrapeg				☐ Hydrapeg			

13	**Multifocal** Patient Age _____ Dominant Eye: ☐ OD ☐ OS	☐ Center Near ☐ Center Distance	Pupil Size (mm)	Spectacle Add	☐ Center Near ☐ Center Distance	Pupil Size (mm)	Spectacle Add

Scan or Photograph and email to DeltaScleral@abboptical.com

Clear Form